

- UK** : Turn on the transmitter and adjust the volume.
D : Empfänger einschalten und Lautstärke einstellen.
NL : Zet ontvanger aan en stel volume in.
FR : Allumez le transmetteur et réglez le volume.
E : Encienda el transmisor y ajuste el volumen.
I : Accendere il trasmettitore e regolare il volume.

1



- UK** : Pull on the underpants.
D : Höschen anziehen.
NL : Broekje aantrekken.
FR : Enfilez le slip.
E : Poner la ropa interior.
I : Togliere le mutandine.

2



- UK** : Attach the transmitter to the underpants.
D : Sender am Höschen befestigen.
NL : Bevestig zender op broekje.
FR : Attachez le transmetteur au slip à l'aide des pressions.
E : Sujetar el transmisor a la ropa interior.
I : Attaccare il trasmettitore alle mutandine.

3



- UK** : Sleep well.
D : Gute Nacht!
NL : Welterusten.
FR : Dormez comme d'habitude.
E : Buenas noches.
I : Dormite bene.

4



- UK** : The transmitter will sound an alarm if urination occurs.
D : Empfänger schlägt beim Wasserlassen Alarm.
NL : Ontvanger geeft alarm bij plassen.
FR : Si vous avez des fuites le transmetteur émettra un signal sonore.
E : El transmisor hará sonar la alarma cuando se vaya a orinar.
I : Il trasmettitore manderà un allarme acustico in caso di urina.

5




- UK** : Get out of bed and go to the toilet.
D : Aufstehen, Wecker ausschalten und zur Toilette gehen.
NL : Opstaan, wekker uit en naar de toilet.
FR : Levez-vous et allez aux toilettes.
E : Salir de la cama e ir al lavabo.
I : Alzatevi e andate in bagno

6



- UK** : Repeat steps 2 and 3. Sleep well.
D : Punkte 2 und 3 wiederholen. Gute Nacht!
NL : Herhaal 2 en 3, weterusten.
FR : Répétez les étapes 2 et 3. Rendornez-vous.
E : Repetir los pasos 2 y 3. Buenas noches.
I : Ripetete le fasi 2 e 3. Dormite bene

7



- UK** : Wake up dry !!!
D : Trocken !!!
NL : Droog !!!
FR : Vous vous réveillerez sec !!!
E : Despertar seco !!!
I : Risveglio asciutto!!!

8

Tips:

- When beginning the therapy, be motivated, relaxed and confident.
- Do not make a game out of the therapy. This will cause stress, and it is counterproductive.
- Never punish your child; your child is not responsible for occasional accidents.
- Do not wake your child at night to go to the toilet.
- Give your child an occasional reward. A small present is a nice positive reinforcement that your child is doing well. Do not stop the therapy.
- Do not start the therapy if important events will soon take place, i.e. moving the family or the birth of a child.
- Leave a small light on; this will keep your child from sleeping too deeply.
- Do not let the temperature of your child's bed get too warm, as this will put your child into a deeper sleep.
- Do not prohibit your child from drinking before going to bed. Drinking less does not decrease the likelihood of bedwetting. However, it is preferable that your child does not drink coffee, tea, or beverages containing carbonation or chocolate before going to bed.
- If your child's bedwetting condition does not improve within four months, speak with your doctor, as there might be a medical problem.

PLEASE READ THIS INFORMATION AND KEEP IT FOR FURTHER USE.

The "RodGer" system is suitable for children over the age of 5 years. Most children stop bedwetting when they are 5. Beyond this age, it might be necessary to help them remain dry at night. Occasionally, a medical disorder may be present, but mostly the bedwetting is due to an involuntary failure of learning; your child is not to blame. Successful treatment with the RodGer system takes from a few weeks to

a few months. Persistence is the key to success. If relapse occurs, another course will usually succeed. Practice this routine before setting the transmitter.

In the early part of treatment, a child who is a heavy sleeper may not wake up. In time, this child will "tune in". Until then, you should wake your child, but make sure that your child stops the alarm and goes to the bathroom.

Use the "RodGer" system every night until there have been no "wet beds" for 14 nights.

To test the system, switch the receiver to the "on" position and keep a metal object at both press-studs of the transmitter.

Keep sufficient ventilation while applying AC power source. To prevent fire or shock hazard, do not expose this product to rain or any type of moisture.

Receiver will release beep sound when the underwear gets wet.

Vibration cushion

For children who are deaf or hard of hearing, a vibration cushion can be ordered

Warranty

We undertake the following guarantee obligations. If within 6 months after receipt of the device by the purchaser manufacturing or material faults are established - at our discretion - these aforementioned faults shall be repaired free of charge. We shall not be liable for more or other than the aforesaid guarantee obligations. The guarantee is only valid if the below-mentioned conditions have been complied with:

- You may only make claim to the aforesaid guarantee if you have the receipt.
- The guarantee obligation expires if the device is used inexpertly or is not used in accordance with the instructions for use or if inexpert repairs are made to the device.
- Batteries and battery lids are not covered by the guarantee.